Leer atentamente:

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Please Read Carefully:

NOTE: THE MEDICAL PROFESSIONAL IS RESPONSIBLE FOR PROPER SURGICAL PROCEDURES AND TECHNIQUES. THE FOLLOWING INSTRUCTIONS ARE FOR INFORMATION ONLY. THE ATTENDING SURGEON MUST EVALUATE THE APPROPRIATENESS OF THE PROCEDURE BASED ON HIS/HER OWN MEDICAL TRAINING AND EXPERIENCE. **THIS PRODUCT IS INTENDED FOR SINGLE USE ONLY.**

Reuse of this device may expose the patient to infection or contamination risks.

Application Instructions

Step 1: IMPORTANT! Cleanse and dry skin on nose and application area.

Step 2: Wipe cleansed area with Isopropyl Alcohol Wipe* and dry thoroughly.

Step 3: Express edema from application area.

Step 4: Apply skin preparation* to application area and DRY THOROUGHLY. VERY IMPORTANT!

Step 5: Apply 1/2 inch wide skin closure tape* starting at the bottom of the nose being careful to overlap each strip to be assured that all the skin in the application area is covered. (See Figure #1)

Step 6: Remove release paper from Dorsum Pad and apply over skin closure tape from Radix to Supratip area. Dorsum Pad adhesive should not contact skin. (See Figure #2)

Step 7: Size the Soft-Form splint to the application area. Splint may be trimmed with scissors if necessary. A properly sized splint should extend laterally to the lateral osteotomy sites, (but not beyond), from the Radix to the Supratip dome junction.

Step 8: Soften nasal splint in 150° -160°F water (approximately 30 seconds).

Step 9: Remove the release paper from the back of the splint and apply adhesive side to nose over the Dorsum Pad and skin closure tape. Mold the Soft-Form™ Splint to achieve proper fit and pressure. **IMPORTANT NOTE:** Do not let splint adhesive contact skin!! (See Figure #3)

Step 10: Splint should be removed in seven to eight days. *Included in the complete kit only.

Removal Instructions

Step 1: Remove by gently spreading blunt bayonet type forceps between the skin closure tape and the patient's skin. Do not use excess force so that the skin tents up.

Step 2: Retaping for three or four days is optional by the physician

Disposal

When disposing of the product, take all steps necessary to avoid risk of injury and infection. Contaminated products must be disposed of as hazardous waste and handled so as to avoid contamination of third parties.

Note

Any serious incident that has occurred in relation to the device should be reported to Invotec International, Inc, MDSS GmbH, and the competent authority of the Member State in which the user and/or patient is established.

Warranty

Invotec International, Inc warrants that the product is free from defects in material and workmanship. Invotec will replace or provide a refund for any product found to be defective so long as the product is returned according to the Returned Goods instructions in the Sales Policy. Invotec shall not be liable for any consequential loss, damage or expense directly or indirectly arising from the use of, or inability to use, this product. **THE FOREGOING WARRANTY IS IN LIEU OF AND EXCLUDES ALL OTHER WARRANTIES, EXPRESSED OR IMPLIED, HOWEVER ARISING, INCLUDING MERCHANTABILITY, FITNESS FOR A PARTICULAR PURPOSE, AGAINST INFRINGEMENT OR OTHERWISE.** Invotec International neither assumes, nor authorizes any person to assume for it, any other additional liability or responsibility with respect to this product.

Soft-Form™ and Custom-Cut are trademarks of Invotec International, Inc.

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Soft-Form™ and Custom-Cut, Invotec International, Inc.'in ticari bir markasıdır.