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Symbols and definitions:

	Manufacturer: Indicates the medical device manufacturer, as defined in EU directives 90/385/EEC, 93/42/EEC, and 98/79/EC.
	Authorized representatives in the European Community. This symbol shall be accompanied by the name and address of the authorized representative in the European Community, adjacent to the symbol.
	Use by date: Indicates the date after which the medical device is not to be used.
	Batch code: Indicates the manufacturer's batch code so that the batch or lot can be identified.
	Catalog Number: Indicates the manufacturer's catalog number so that the medical device can be identified.
	Do not use if package is damaged: Indicates a medical device that should not be used if the package has been damaged or opened.
	Do not re-use: Indicates a medical device that is intended for one use, or for use on a single patient during a single procedure.
	Consult instructions for use: Indicates the need for the user to consult instructions for use.
	Caution: Indicates the need for the user to consult the instructions for use for important cautionary information such as warnings and precautions that cannot, for a variety of reasons, be presented on the medical device itself.
	Medical Device Indicated

Please Read Carefully:

NOTE: THE MEDICAL PROFESSIONAL IS RESPONSIBLE FOR PROPER SURGICAL PROCEDURES AND TECHNIQUES. THE FOLLOWING INSTRUCTIONS ARE FOR INFORMATION ONLY. THE ATTENDING SURGEON MUST EVALUATE THE APPROPRIATENESS OF THE PROCEDURE BASED ON HIS/HER OWN MEDICAL TRAINING AND EXPERIENCE.

**This product is intended for single use only.** Reuse of this device may expose the patient to infection or contamination risks.

**Application Instructions**

Step 1: IMPORTANT...Cleanse and dry skin on nose and application area.

Step 2: Wipe cleansed area with Isopropyl Alcohol Wipe\* and dry thoroughly.

Step 3: Express edema from application area.

Step 4: Apply skin preparation\* to application area and DRY THOROUGHLY... VERY IMPORTANT!

Step 5: Apply 1/2 inch wide skin closure tape\* starting at the bottom of the nose being careful to overlap each strip to be assured that all the skin in the application area is covered. (See Figure #1)

Step 6: Remove release paper from Dorsum Pad and apply over skin closure tape from Radix to Supratip area. Dorsum Pad adhesive should not contact skin. (See Figure #2)

Step 7: Size the Soft\_Form splint to the application area. Splint may be trimmed with scissors if necessary. A properly sized splint should extend laterally to the lateral osteotomy sites, (but not beyond), from the Radix to the Supratip dome junction.

Step 8: For easier removal, pull up a small corner of the backing before putting the splint into the hot water. (150°-160°F/65°-75°C and 30 seconds or when splint is transparent)

Step 9: Remove the release paper from the back of the splint and apply adhesive side to nose over the Dorsum Pad and skin closure tape. Mold the Soft-Form Splint to achieve proper fit and pressure. IMPORTANT NOTE: Do not let splint adhesive contact skin!! (See Figure #3)

Step 10: Splint should be removed in seven to eight days.

\*Included in the complete kit only.

**Removal Instructions**

Step 1: Remove by gently spreading blunt bayonet type forceps between the skin closure tape and the patient's skin. Do not use excess force so that the skin tents up.

Step 2: Retaping for three or four days is optional by the physician

**Disposal**

When disposing the product, take all steps necessary to avoid risk of injury and infection. Contaminated products must be disposed of as hazardous waste and handled so as to avoid contamination of third parties.

**Note**

Any serious incident that has occurred in relation to the device should be reported to Invotec International, Inc, MDSS GmbH, and the competent authority of the Member State in which the user and/or patient is established.

**Warranty**

Invotec International, Inc. warrants that the product is free from defects in material and workmanship. Invotec will replace or provide a refund for any product found to be defective so long as the product is returned according to the Returned Goods instructions in the Sales Policy. Invotec shall not be liable for any consequential loss, damage or expense directly or indirectly arising from the use of, or inability to use, this product. THE FOREGOING WARRANTY IS IN LIEU OF AND EXCLUDES ALL OTHER WARRANTIES, EXPRESSED OR IMPLIED, HOWEVER ARISING, INCLUDING MERCHANTABILITY, FITNESS FOR A PARTICULAR PURPOSE, AGAINST INFRINGEMENT OR OTHERWISE. Invotec International neither assumes, nor authorizes any person to assume for it, any other additional liability or responsibility with respect to this product. Soft-Form and Custom-Cut are trademarks of Invotec International, Inc.

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Soft-Form ve Custom-Cut Invotec International, Inc. ticari markasıdır.