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6833 Phillips Industrial Blvd. • Jacksonville, Florida 32256 U.S.A.

Tel: 800-998-8580 / 904-880-1229 • Fax: 904-886-9517 • www.invotec.net

PJ SMITH MEDICAL PRODUCTS MARKETING

18 Yeates Close • Thame OX93AR, UK

INVOTEC INTERNATIONAL, INC.

6833 Phillips Industrial Blvd. • Jacksonville, Florida 32256 U.S.A.

Tel: 800-998-8580 / 904-880-1229 • Fax: 904-886-9517 • www.invotec.net

PJ SMITH MEDICAL PRODUCTS MARKETING

18 Yeates Close • Thame OX93AR, UK

Please Read Carefully:

NOTE: THE MEDICAL PROFESSIONAL IS RESPONSIBLE FOR PROPER SURGICAL PROCEDURES AND TECHNIQUES. THE FOLLOWING INSTRUCTIONS ARE FOR INFORMATION ONLY. THE ATTENDING SURGEON MUST EVALUATE THE APPROPRIATENESS OF THE PROCEDURE BASED ON HIS/HER OWN MEDICAL TRAINING AND EXPERIENCE.

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The Quick-Pack Nasal Catheter was designed as a packing following septoplasty, rhinoplasty, or other intra-nasal surgical procedures. It may also be used as nasal packing to control anterior nasal hemorrhage.

#### How To Use The Quick-Pack Catheter.

Step 1: Test inflate the cuff before introducing the catheter. The Inflation valve is marked to indicate its volume.

Step 2: Lubricate the catheter with sterile water or non-petroleum based lubricant.

Step 3: Remove blood clots from nose.

Step 4: Insert the catheter so that the cuff rest in the anterior chamber.

Step 5: Inflate the cuff with either sterile water or sterile saline by inserting a syringe in the valve marked 30 CC. A 20 CC syringe is provided. NOTE: Inflation volumes of 10 to 25 CC's are usually adequate.

Step 6: Check cuff pressure frequently and relieve cuff pressures when indicated. NOTE: A small well lubricated suction catheter may be inserted through the airway to remove secretions and/or blood clots.

#### Removal

Insert syringe firmly into valve and evacuate cuff completely. After verification of complete cuff deflation, gently withdraw catheter.

#### To Sterilize

If re-sterilization is necessary, STEAM AUTOCLAVING is the recommended method. The Ultra-Stat catheter should be washed with mild soap, thoroughly rinsed, wrapped in a lint free wrapper and steam autoclaved at 121° C ( 250°F) for 30 minutes.

THIS PRODUCT IS FOR SINGLE USE ONLY.

*Caution: U.S. Federal Law restricts this device to sale by or on the order of a physician.*

#### Warranty

Invotec International, Inc. warrants that the product is free from defects in material and workmanship. Invotec will replace or provide a refund for any product found to be defective so long as the product is returned according to the Returned Goods instructions in the Sales Policy. Invotec shall not be liable for any consequential loss, damage or expense directly or indirectly arising from the use of, or inability to use, this product. THE FOREGOING WARRANTY IS IN LIEU OF AND EXCLUDES ALL OTHER WARRANTIES, EXPRESSED OR IMPLIED, HOWEVER ARISING, INCLUDING MERCHANTABILITY, FITNESS FOR A PARTICULAR PURPOSE, AGAINST INFRINGEMENT OR OTHERWISE. Invotec International neither assumes, nor authorizes any person to assume for it, any other additional liability or responsibility with respect to this product.

*Quick-Pack is a trademark of Invotec International, Inc.*

INVOTEC INTERNATIONAL, INC.

6833 Phillips Industrial Blvd. • Jacksonville, Florida 32256 U.S.A.

Tel: 800-998-8580 / 904-880-1229 • Fax: 904-886-9517 • www.invotec.net

PJ SMITH MEDICAL PRODUCTS MARKETING

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*Quick-Pack Invotec International Inc'nin bir ticari markasıdır.*