BOTÓN SEPTAL

Leer atentamente:

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SEPTAL BUTTON

Please Read Carefully:

NOTE: THE MEDICAL PROFESSIONAL IS RESPONSIBLE FOR PROPER SURGICAL PROCEDURES AND TECHNIQUES. THE FOLLOWING INSTRUCTIONS ARE FOR INFORMATION ONLY. THE ATTENDING SURGEON MUST EVALUATE THE APPROPRIATENESS OF THE PROCEDURE BASED ON HIS/HER OWN MEDICAL TRAINING AND EXPERIENCE.

THIS PRODUCT IS INTENDED FOR SINGLE USE ONLY.

Reuse of this device may expose the patient to infection or contamination risks.

Indications

When a patient is a poor candidate for surgical correction of a septal perforation, and the nature of the perforation, that is size and location, and septal anatomy are appropriate, the attending surgeon may consider closure with a septal button. Generally, anterior septal perforations that are 2cm. or less in maximum diameter are suitable for closure with a septal button. This product is for single use only. Reuse of this device may expose the patient to infection or contamination risks. Once the seal of the sterile packaging has been torn open, its contents will not be taken back by the manufacturer.

Contraindications

Septal deformities are relative contraindications for septal button placement, as the flange on the convex side of the deformity may increase nasal obstruction. Relative contraindications may also include active local nasal infection as the septal button may prevent adequate treatment.

How to use the Septal Button

Step 1: After adequate topicalization of the nose, 1% lidocaine with 1:100,000 epinephrine is infiltrated about the margins of the perforation. Light sedation may be required in some cases.

Step 2: Custom trim the flanges of the button to permit 3-4mm. of overlap with respect to the margins of the perforation.

Step 3: Choose one of the following two methods of insertion:

METHOD 1:

Collapse one flange and introduce the button into the perforation, then gently pull it back through from the opposite end.

METHOD 2:

Sew the edges of one flange together in a rosette manner with a suture, tie suture and leave long. Pass the long end of the suture through the perforation and gently guide the button in place using a hemostat from the back side to apply gentle pressure. Cut and remove the suture allowing the sewn flange to open flat.

NOTE: CARE MUST BE TAKEN SO THAT THE FLANGES OF THE BUTTON DO NOT ABUTT THE POSITION OF THE SEPTUM WITH THE UPPER LATERAL CARTILAGE, AS THIS IS A PARTICULARLY ANNOYING SENSATION FOR THE PATIENT.

Removal

Remove button by carefully cutting the center post with a scalpel or microscissors and removing each piece of the button through separate nostrils.

Caution: U.S. Federal law restricts this device to sale by or on the order of a physician.

Disposal: When disposing the product, take all steps necessary to avoid risk of injury or infection. Contaminated products must be disposed of as hazardous waste and handled to avoid contamination of third parties.

Warranty

Invotec International, Inc. warrants that the prAoduct is free from defects in material and workmanship. Invotec will replace or provide a refund for any product found to be defective so long as the product is returned according to the Returned Goods instructions in the Sales Policy. Invotec shall not be liable for any consequential loss, damage or expense directly or indirectly arising from the use of, or inability to use, this product. THE FOREGOING WARRANTY IS IN LIEU OF AND EXCLUDES ALL OTHER WARRANTIES, EXPRESSED OR IMPLIED, HOWEVER ARISING, INCLUDING MERCHANTABILITY, FITNESS FOR A PARTICULAR PURPOSE, AGAINST INFRINGEMENT OR OTHERWISE. Invotec International neither assumes, nor authorizes any person to assume for it, any other additional liability or responsibility with respect to this product.

Invotec International neither assumes, nor authorizes any person to assume for it, any other additional liability or responsibility with respect to this product.

This product is for single use only.

Reuse of this device may expose patient to infection or contamination risks. once the seal of the sterile packaging has been torn open its contents will not be taken back by the manufacturer. Any serious incident that has occurred in the relation to the device should be reported to Invotec international, Inc. MDSS GmbH, and the competent authority of the Member State in which the user and/or patient is established.

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