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Symbols and definitions:


	Manufacturer: Indicates the medical device manufacturer, as defined in EU directives 90/385/EEC, 93/42/EEC, and 98/79/EC.
	Authorized representatives in the European Community. This symbol shall be accompanied by the name and address of the authorized representative in the European Community, adjacent to the symbol.
	Use by date: Indicates the date after which the medical device is not to be used.
	Batch code: Indicates the manufacturer's batch code so that the batch or lot can be identified.
	Catalog Number: Indicates the manufacturer's catalog number so that the medical device can be identified.
	Sterilized using ethylene oxide: Indicates a medical device that has been sterilized using ethylene oxide.
	Do not resterilize: Indicates a medical device that is not to be resterilized.
	Do not use if package is damaged: Indicates a medical device that should not be used if the package has been damaged or opened.
	Do not re-use: Indicates a medical device that is intended for one use, or for use on a single patient during a single procedure.
	Consult instructions for use: Indicates the need for the user to consult instructions for use.
	Caution: Indicates the need for the user to consult the instructions for use for important cautionary information such as warnings and precautions that cannot, for a variety of reasons, be presented on the medical device itself.

REV 112618

TRACHEAL T-TUBE

PATIENT INFORMATION

POST-OPERATIVE CARE AND DISCHARGE INSTRUCTIONS

Application

Tracheal T-Tubes are designed to keep the cervical and thoracic (upper and middle) tracheal airway patent and to assist ventilation through the larynx. They can allow normal phonation and humidification and may provide support in the management of acute tracheal injury.

Maintenance and Hygiene

For the patient with a Tracheal T-Tube, post-op care and proper instructions are vital to the successful application of the device.

The T-Tube must be kept patent and clean. At home, the patient may benefit from a room humidifier, especially in the bedroom. To avoid losing the outer (capping) plug, the patient may pass a thread through the plug, make a loop and then secure the thread with a pin or tape. Patient should always have an extra Tracheostomy tube available.

Post-op Instructions, Week 1

1. Keep small tracheostomy tube at patient's bedside for physician or emergency use.
2. Plug external branch of T-Tube when patient is recovering from anesthesia. Unplug for any signs of respiratory distress.
3. Constant humidity is essential to keep secretions moist and to help in mobilizing secretions. It may be necessary to remove the plug and apply humidity directly over the external opening. Use a face mask if tube is plugged. Use a tracheostomy mask if tube is unplugged.
4. Instill saline through external branch to prevent and dislodge crusting within the T-Tube, three times daily.
5. Clean skin around external portion of T-Tube with cottontipped applicator and half-strength hydrogen peroxide, three times daily. Apply Providone iodine ointment after cleaning.
6. Clean inside of external branch of T-Tube with cotton-tipped applicator dipped in half-strength hydrogen peroxide.

Post-op Instructions, after Week 1 and for Discharge

1. Continue plugging external branch of T-Tube as directed. Normal respiration should keep T-Tube and secretions moist.
2. Instill T-Tube with saline two times daily.
3. Maintain skin care twice daily.
4. Clean external branch of T-Tube twice daily.

Procedure for Instillation with saline

1. Unplug T-Tube, instill 1/2 to 1 cc of saline. Never use syringe with needle attached when instilling.
2. Plug T-Tube and have patient cough to clear secretions. Patient may require suctioning at this time.
3. Repeat step 1 as needed with up to 1 to 3cc of normal saline. Repeat step 2

Procedure for Suctioning.

1. Remove plug and suction:
 - Upwards: the suction catheter is directed up the T-Tube by bending the external branch down
 - Downwards: the suction catheter is directed down the T-Tube by bending the external branch upIt is not necessary to insert the suction catheter much beyond the end of the T-Tube. This could produce irritation of the tracheal lining and coughing. Deep suction only when the patient's condition warrants.

Homecare Method for removing crusts

1. Place mirror and light on countertop.
2. Set up suction machine so that tip is readily available.
3. While seated, adjust mirror and light to provide an adequate view of the inside external branch of the T-Tube.
4. Using small alligator forceps, carefully manipulate pincers to grasp crusts or plugs and withdraw; use suction to remove smaller pieces.
5. Repeat until the inside of the tube is completely clear Note: Some patients experience excessive crust formation and the formation of mucous plugs. Use of an Ultrasonic Humidifier, a Nebulizer Machine or an Atomizer Machine can improve the situation. Please check with your physician or other Healthcare professional for any other specific needs.

WARRANTY

Invotec International, Inc. warrants that the product is free from defects in material and workmanship. Invotec will replace or provide a refund for any product found to be defective so long as the product is returned according to the Returned Goods instructions in the Sales Policy. Invotec shall not be liable for any consequential loss, damage or expense directly or indirectly arising from the use of, or inability to use, this product. THE FOREGOING WARRANTY IS IN LIEU OF AND EXCLUDES ALL OTHER WARRANTIES, EXPRESSED OR IMPLIED, HOWEVER ARISING, INCLUDING MERCHANTABILITY, FITNESS FOR A PARTICULAR PURPOSE, AGAINST INFRINGEMENT OR OTHERWISE. Invotec International neither assumes, nor authorizes any person to assume for it, any other additional liability or responsibility with respect to this product.

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