ULTRA STAT™
CATÉTER PARA EPISTAXIS

Leer atentamente:

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ULTRA-STAT™
EPISTAXIS CATHETER

Please Read Carefully:

NOTE: THE MEDICAL PROFESSIONAL IS RESPONSIBLE FOR PROPER SURGICAL PROCEDURES AND TECHNIQUES. THE FOLLOWING INSTRUCTIONS ARE FOR INFORMATION ONLY. THE ATTENDING SURGEON MUST EVALUATE THE APPROPRIATENESS OF THE PROCEDURE BASED ON HIS/HER OWN MEDICAL TRAINING AND EXPERIENCE.

Indications

The Ultra-Stat Catheter was designed as a nasal packing for controlling nasal hemorrhage. Ultra-Stat may also be used for nasal packing following septoplasty, rhinoplasty, or other intra-nasal surgical procedures.

Warnings

- Do not over inflate the cuff. Overinflation may partially obstruct the airway causing restricted or labored breathing, or may result in tissue necrosios due to excessinve pressure.

- Do not exceed the volumes recommended for filling the cuff. Exceeding these volumes may cause the cuffs to rupture and cause trauma to the surrounding intranasal tissues or dislodgement of the nasal catheter.

How to use the Ultra-Stat Balloon Catheter

- Test inflate the cuff before introducing the catheter. Each inflation valve is marked to indicate its corresponding cuff by volume.
- Lubricate the catheter with sterile water or non-petroleum based lubricant.
- Remove blood clots from nose.
- Insert the catheter so that the bevelled tip rest in the nasopharynx.
- Inflate the nasopharyngeal cuff with either sterile water or sterile saline by inserting the syringe in the valve marked 10 CC. A 20 CC syringe is provided.
- Remove the syringe and gently pull the catheter to seat it.
- Inflate the anterior balloon in the same manner using up to 30 CC’s of sterile solution. NOTE: Inflation of 10 to 25 CC’s are usually adequate.
- Check both cuff pressures frequently and relieve cuff pressures when indicated. NOTE: A small, well lubricated suction catheter may be inserted through the airway to remove secretions and/or blood clots.

Removal of Catheter

Insert syringe firmly into each valve, in turn, and evacuate each cuff completely. After verification of complete cuff deflation, gently withdraw catheter.

Disposal: When disposing the product, take all steps necessary to avoid risk of injury or infection. Contaminated products must be disposed of as hazardous waste and handled to avoid contamination of third parties.

General Handling

Any contamination of the product must be avoided.

The product must be kept in its sealed protective packaging.

Examine packaging for damage before opening sdnce damaged packaging may impair the sterility of the product.

Do not open the protective until just before implantation of the product. In addition, the product must be visually checked for damage.

This product is intended for single use only. Reuse of this device may expose patient to infection or contamination risks. once the seal of the sterile packaging has been torn open its contents will not be taken back by the manufacturer.

THIS PRODUCT IS FOR SINGLE USE ONLY.

Any serious incident that has occurred in the relation to the device should be reported to Invotec international, Inc. MDSS GmbH, and the competent authority of the Member State in which the user and/or patient is established.

Caution: U.S. Federal law restricts this device to sale by or on the order of a physician.

Warranty

Invotec International, Inc. warrants that the prAoduct is free from defects in material and workmanship. Invotec will replace or provide a refund for any product found to be defective so long as the product is returned according to the Returned Goods instructions in the Sales Policy. Invotec shall not be liable for any consequential loss, damage or expense directly or indirectly arising from the use of, or inability to use, this product. THE FOREGOING WARRANTY IS IN LIEU OF AND EXCLUDES ALL OTHER WARRANTIES, EXPRESSED OR IMPLIED, HOWEVER ARISING, INCLUDING MERCHANTABILITY, FITNESS FOR A PARTICULAR PURPOSE, AGAINST INFRINGEMENT OR OTHERWISE. Invotec International neither assumes, nor authorizes any person to assume for it, any other additional liability or responsibility with respect to this product.

Ultra-Stat is a trademark of Invotec International, Inc.

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Ultra-Stat, Invotec International, Inc.'in ticari bir markasıdır.
6833 Phillips Industrial Blvd. • Jacksonville, Florida 32256 U.S.A.