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**Product description**

Ventilation tube for ventilation of the middle ear. The tube should be implanted via a small incision in the eardrum (myringotomy). The tube has a central lumen of varying diameter. The actual size to use should be decided by the surgeon based on the respective indication. The tube facilitates the removal of secretions (middle ear effusion) and relieving the auditory tube (Eustachian tube).

**Indications**

- Tympanic effusion (secretory otitis media, glue ear)
- Chronic ventilation problems

**Contraindications**

Allergic reactions to the material used.

**Side Effects**

1. In rare cases otitis media is possible after myringotomy; in such cases the patient should be treated with an antibiotic.
2. Insertion of an ontological ventilation tube very rarely results in a potentially permanent eardrum perforation in which in turn can be closed through tympanoplastic surgery, if necessary.
3. The tube may also fall into the middle ear. In this case the tube must be surgically removed.
4. Usually the tube will be rejected spontaneously after a prolonged period, when it simply falls out onto the outer ear. If this does not occur then the tube should be surgically removed, though this is not normally required.
5. If strong viscid secretions are produced (secretory otitis media) the tube may become clogged. If it is not possible to free the tube from these clogging secretions then a new vent tube must be inserted.
6. Premature rejection of the ventilation tube may occur. If this occurs, new ontological ventilation tube may have to be inserted after the reasons for the rejection have been identified.
7. There is the possibility of scarring of formation of hyaline structure on eardrum level.
8. Improper insertion may cause luxation of the ossicular chain.
9. The existence of anatomical anomalies may lead to injury to the bulbous vena jugularis.

**Precautions**

While the tube remains inserted in the eardrum the patient should not allow water to flow through the tube into the middle ear, since this may result in a secondary infection. This can occur as a result of bathing, washing the hair, swimming, etc. To avoid such infection the use of tight fitting ear plugs is recommended. The patient should also avoid any air pressure variations while

the tube remains in situ. All ventilation tubes from Invotec International are MRI compatible to 4.5T.

**General handling**

Any contamination of the product must be avoided. The product must be kept in its sealed protective packaging. Examine packaging for damage before opening since damaged packaging may impair the sterility of the product. Do not open the protective packaging until just before implantation of the product. In addition, the product must be visually checked for damage. This product is intended for single use only. Once the seal of the sterile packaging has been torn open its contents will not be taken back by the manufacturer.

**Packaging and sterility**

The packaging of the product consists of • Sterile packaging (primary packaging) The packaging meets the provisions of the European Standards. Intact packaging protects the product from environmental influences and ensures sterile storage.

**Handling of sterile packaging**

Please ensure that the relevant aseptic instructions are complied with when removing the product from the packaging.

**Re-sterilization**

Re-sterilization after expiry of the use-by date as well as re-sterilization of unsterile uncontaminated products is prohibited. Manufacturer and distributor will assume no liability for products that have been re-sterilized by the user.

**Storage**

Store in a dry environment at room temperature. After expiry of the use-by date the product may no longer be used.

**Warranty**

Invotec International, Inc. warrants that the product is free from defects in material and workmanship. Invotec will replace or provide a refund for any product found to be defective so long as the product is returned according to the Returned Goods instructions in the Sales Policy. Invotec shall not be liable for any consequential loss, damage or expense directly or indirectly arising from the use of, or inability to use, this product. THE FOREGOING WARRANTY IS IN LIEU OF AND EXCLUDES ALL OTHER WARRANTIES, EXPRESSED OR IMPLIED, HOWEVER ARISING, INCLUDING MERCHANTABILITY, FITNESS FOR A PARTICULAR PURPOSE, AGAINST INFRINGEMENT OR OTHERWISE. Invotec International neither assumes, nor authorizes any person to assume for it, any other additional liability or responsibility with respect to this product.

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\*Hydromer est une marque de fabrication de Hydromer, Inc.