Leer atentamente:

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Schiffgraben 41 • 30175 Hannover, Germany

Please Read Carefully:

NOTE: THE MEDICAL PROFESSIONAL IS RESPONSIBLE FOR PROPER SURGICAL PROCEDURES AND TECHNIQUES. THE FOLLOWING INSTRUCTIONS ARE FOR INFORMATION ONLY. THE ATTENDING SURGEON MUST EVALUATE THE APPROPRIATENESS OF THE PROCEDURE BASED ON HIS/HER OWN MEDICAL TRAINING AND EXPERIENCE.

Indications

The Post-Stat nasal balloon catheter was designed as a posterior nasal packing for controlling posterior nasal hemorrhage and as an inter-operative packing to aid in the prevention of fluid aspiration. The Post-Stat catheter allows for both suction and irrigation.

Warnings

- Do not over inflate the cuff. Overinflation may partially obstruct the airway causing restricted or labored breathing, or may result in tissue necrosis due to excessive pressure.
- Do not exceed the volumes recommended for filling the cuff. Exceeding these volumes may cause the cuffs to rupture and cause trauma to the surrounding intranasal tissues or dislodgement of the nasal catheter.

How To Use The Post-Stat Balloon Catheter

- 1: Test the cuff before introducing the catheter by inflating with 5 to 10 CC’s of sterile saline solution or sterile water.
- 2: Lubricate the catheter with sterile water or non-petroleum based lubricant.
- 3: If blood clots are present remove the clots from the nose.
- 4: Insert the catheter so that the beveled tip rest in the nasopharynx.
- 5: Withdraw the stylet.
- 6: Inflate the nasopharyngeal cuff with either sterile water or sterile saline by inserting a syringe in the inflation valve. A 20 CC syringe is provided. NOTE: Inflation volumes of 4 to 8 CC’s are usually adequate.
- 7: Remove the syringe and gently pull the catheter to seat it.
- 8: If suction and/or irrigation is to be used set up proper connectors.
- 9: Check cuff pressure frequently and relieve cuff pressure when indicated.

Removal

Insert syringe firmly into valve, and evacuate cuff completely. After verification of complete cuff deflation, gently withdraw catheter. **Disposal:** When disposing the product, take all steps necessary to avoid risk of injury or infection. Contaminated products must be disposed of as hazardous waste and handled to avoid contamination of third parties.

THIS PRODUCT IS FOR SINGLE USE ONLY.

Reuse of this device may expose patient to infection or contamination risks. once the seal of the sterile packaging has been torn open its contents will not be taken back by the manufacturer. Any serious incident that has occurred in the relation to the device should be reported to Invotec international, Inc. MDSS GmbH, and the competent authority of the Member State in which the user and/or patient is established.

Caution: U.S. Federal Law restricts this device to sale by or on the order of a physician.

Warranty

Invotec International, Inc. warrants that the product is free from defects in material and workmanship. Invotec will replace or provide a refund for any product found to be defective so long as the product is returned according to the Returned Goods instructions in the Sales Policy. Invotec shall not be liable for any consequential oss, damage or expense directly or indirectly arising from the use of, or inability to use, this product. THE FOREGOING WARRANTY IS IN LIEU OF AND EXCLUDES ALL OTHER WARRANTIES, EXPRESSED OR IMPLIED, HOWEVER ARISING, INCLUDING MERCHANTABILITY, FITNESS FOR A PARTICULAR PURPOSE, AGAINST INFRINGEMENT OR OTHERWISE. Invotec International neither assumes, nor authorizes any person to assume for it, any other additional liability or responsibility with respect to this product.

Post-Stat is a trademark of Invotec International, Inc.

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Ultra-Stat, Invotec International, Inc.'in ticari bir markasıdır.
6833 Phillips Industrial Blvd. • Jacksonville, Florida 32256 U.S.A.