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Schiffgraben 41 • 30175 Hannover, Germany

0197

These instructions are for professional reference prior to surgery and should be read carefully.

#### Description

Invotec Tracheal T-Tubes are designed for the surgical management of tracheal and subglottic stenosis and the reconstruction of cervical and thoracic trachea. The T-Tubes, made of medical grade silicone, are designed to maintain patency of the tracheal airway and to provide ventilation through the larynx. They can allow normal humidification and phonation and may provide support in the management of acute tracheal injury.

#### Indications

The T-Tube may be used as both a tracheal stent and a tracheostomy tube and may be used under the following circumstances:

- In tracheal stenosis, when the cervical or thoracic airway cannot be repaired, or as a substitute for the cervical trachea when it cannot be reconstituted or reconstructed;
- As a support following reconstructed cervical trachea;
- Prior to reconstruction to maintain an adequate airway while waiting for inflammation to subside;
- As a palliative measure for patients with unresectable carcinoma of the trachea;
- When there is not enough trachea left to repair because of prior surgery;
- With segmental resection and anastomosis;
- When the patient is not a candidate for reconstructive tracheal surgery.

#### Contraindications

The T-Tube cannot be used to prevent aspiration or used as a conduit for positive ventilation through the external branch. It should be predetermined if the proximal limb will rest on the vocal cords with additional risks discussed by health care personnel.

#### Caution:

- This product is for single patient use only
- Federal Law restricts the sale of this device by or on the order of a physician
- This product is shipped NON-STERILE

#### Sterilization

**Step 1.** Remove product from package; using clean gloves, rinse T-Tube with clean water. Pat dry with lint free towel.

**Step 2.** Place product in sterilization package and record proper information on sterilization label.

**Step 3.** Sterilize using Steam Autoclave: Steam Sterilization: flash sterilize using steam for 3 minutes at 270 degrees F at 30 psi.

#### Handling Techniques

1. The Tracheal T-Tube should remain plugged or have a speaking valve attached when possible to help keep the trachea moist and allow normal phonation and respiration.

2. Before cleaning, remove the plug. To avoid misplacement of the plug, set aside in a clean place.

3. The Tracheal T-Tube may be suctioned 2 to 3 times per day or as directed by the physician. Tilt external portion of tube upward or downward to direct suction tip to upper or lower ends. Insert suction catheter slowly and carefully while removing mucus and secretions.

4. Clean inside and outside of external branch with cotton-tipped applicator dipped in hydrogen peroxide. Using another applicator dipped in saline solution, continue to clean until no hydrogen peroxide is left on applicator and the inside of the Tracheal T-Tube external branch is clear.

5. Twice daily, apply Betadine® solution or antibiotic ointment to clean skin around external portion of tube.

6. During the winter months, or when in prolonged indoor living conditions with automatic climate control, to prevent occurrence persistent crusting occurs inside the tube, constant humidity is recommended.

- The Intra-operative placement during an open resection-reanastomosis procedure is achieved by careful positioning of the T-Tube to support the freshly reconstructed site. The external branch (horizontal arm) is ideally placed away from the anastomosis to minimize the risk of healing delay.
- Endoscopic placement of the T-Tube can be difficult and tedious. The correct T-Tube diameter should be determined so as to allow for little or no horizontal movement. In cases of subglottic stenosis the T-Tube is accurately tailored so that it does not extend above the free border of the vocal cords. The T-Tube is preferably trimmed to suit the precise dimensional requirements.

1. The placement of the T-Tube is preceded by dilation and calibration.

2. On the basis of the calibrated diameter, a rubber dilator, flexible catheter or umbilical tape may be introduced through the tracheostoma and passed retrograde through the larynx, grasped using a laryngoscope and drawn out through the mouth.

3. The T-Tube is correctly oriented and placed either directly into the stoma and pulled retrograde with the external limb secured or attached extra-orally and transported prograde into the trachea until the external limb can be extracted through the stoma, allowing appropriate positioning of the superior an inferior T-Tube shaft limbs within the trachea.


4. When completed, the T-Tube inferior limb is distal to the stoma and the superior limb extends proximally, with the external limb protruding from the stoma.

5. The dilators, catheters or umbilical tape are disengaged and the T-Tube position checked for accurate level placement.

- 6. An endoscopic examination can verify that the T-Tube is correctly positioned without buckling and the patient is allowed to awaken. Improper length of the T-Tube is corrected as soon as recognized.
- The use of the flexible dilators in a retrograde passage technique allows the stenotic segment(s) to be atrautomatically dilated and the T-Tube drawn into the stenotic segment for support. The horizontal arm is plugged and the patient's airway status evaluated. When a larger size T-Tube is desirable, the endoscopic procedure is repeated after several weeks and larger caliber T-Tube is inserted. When the final lumen size has been achieved, the patient is evaluated for final T-Tube removal which is performed under anesthesia to permit endoscopic evaluation.

- A properly placed T-Tube may be left in-situ up to 29 days to allow firm healing and stabilization of the repair.
- Removal of the T-Tube is best performed under anesthesia so that the airway patency can be evaluated endoscopically, undue patient discomfort is avoided, and granulations, sutures and other impediments may be removed.

Before removing the T-Tube, suction completely and confirm that there is no evidence of excess deposits on the interior of the T-Tube.

TRACHEAL T-TUBE  
OPERATING ROOM PERSONNEL

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6833 Phillips Industrial Blvd. • Jacksonville, Florida 32256 U.S.A.
MDSS GmbH