



# Brightstar Loupe

## Binocular Loupe Headlight Attachment

The high quality, coated optics ensures high resolution and color authenticity. The angle of convergence of the telescopic spectacles can be aligned independently of the inter-pupil distance and allows the optimal adjustment. The flip-up magnifying unit enables spatial orientation during work pauses.



### Technical Data

<b>Item Code</b>	<b>34-20200N</b>	<b>34-20205N</b>
<b>Illumination</b>	F/O	F/O
<b>Magnification</b>	1,8x	2,5x
<b>Working Distance (mm)</b>	370	400
<b>Field of View (mm)</b>	90	80